

K-9 BSD-2 Reward Use & Info

The K-9 BSD-2 uses include balls with 2.75 diameters or smaller, large Kong 2.75" x 4.5" or Terry Towel 13"x 17" folded twice (lengthwise) and rolled.



Towel, Kong, or ball with string is recommended to avoid risk of dog swallowing or chasing into hazardous setting.

K-9 BSD Odor Imprinting Ball



Rubber balls, tennis balls & Kong's have an abundant odor of their own. Our **Odor Imprinting Balls** have minimal odor of their own, allowing the odor of the training aid you put inside to be the most abundant odor source. These balls have holes evenly distributed in them to allow the odor to be easily detected by your K-9. These balls are similar in size to a tennis ball, made extra thick from a durable poly plastic. There is one elongated slot on the ball to allow you to place your training aid (marijuana, explosive scent, cadaver scent, etc.) inside. Place your training aid into small cotton bag (*part of BSD-2 Kit*) and insert into the ball. You can also use odor impregnated cotton balls and place them inside the ball. There is a 12" rope attached to the ball so you can combine the fight (tug action) with the odor reward.

Training protocol: Dog is conditioned to search cue such as "find It" by hiding in cover, debris pile, vehicle etc. dog is then sent off lead with timely praise and play for locating. Only when dog is conditioned to the cue and showing consistent hunt drive to odor source (odor recognition) is the alert behavior trained.